



# ELECTRIC MEDICAL EXAMINATION COUCH

P.C. No.: 1240M

With Gas Spring



#### STANDARD FEATURES

- 2 sectional Top made 16 Gauge M.S. ERW pipes.
- 2 sectional padded Top with water Resistant vinyl fabric.
- Back UP-DOWN on Hydraulic Gas spring.
- Electrically hi-low system on actuator with wired hand set, provided with manual CPR on crank.
- Base trolley made from 16 gauge.
- Base trolley mounted on 125 mm dia plastic body central locking castors.
- One swival B.P. tray.
- One stainless steel I.V. pole.

#### TECHNICAL SPECIFICATIONS\*

Approximate

Overall length	203 cms (80")
Overall width	63 cms (25")
Height min.	63cm (25")
Height max.	81 cms (32")
Sleep deck area	203 x 63 cms (80" x 25")

### **CONSULTING & OPD EQUIPMENT**

# MEDICAL EXAMINATION COUCH

P.C. No.: 1240 Back raising without Gas Spring

P.C. No.: 1240G Back raising with Gas Spring

# EXAMINATION TABLE

P.C. No.: 1241 Made of Heavy Gauge M.S. Tubes

P.C. No.: 1241S Made of 304 Stainless Steel





#### STANDARD FEATURES

- Made of 18G CRC sheet construction.
- 3 pull drawers, center drawer with lock featurs.
- 2 No's lockable sliding doors movement on stainless steel channel.
- One inner shelf full length for storage .
- One swivel BP instrument tray.
- One stainless steel I.V. pole.
- One special built in step stool recessable.
- Heavy duty padded top with water resistant vinyl fabric.
- Epoxy, powder coated finished.

- Heavy padded 2 section top with water resistant vinyl fabric.
- Back section operated with gas spring.
- One paper roller holder.
- Knock down construction.
- Epoxy, powder coated finished.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	200 cms (80")
Overall width	63.5 cms (25")
Overall height	84 cms (33")

TECHNICAL SPECIFICATIO	NS* Approximate
Overall length	195.5 cms (77")
Overall width	62 cms (24 ½")
Overall height	75 cms (29 ½")
Overall dimension	195 cms X 62 cms X 75 cms
	(77" X 24½" X 29½")



## REVOLVING HEIGHT ADJUSTABLE STOOL

P. C. No.: 1270

S.S. top with gas spring without castors



#### **REVOLVING STOOL**

P. C. No.: 1270A

Height adjustable with screw



#### DOCTORS REVOLVING EXAMINATION HEIGHT ADJUSTABLE STOOL

P. C. No.: 1271



#### STANDARD FEATURES

- Height adjustable with gas spring.
- Plastic non rusting base with five rubber shoes.
- Stainless steel top.

TECHNICAL SPECIFICATION	NS* Approximate
Minimum height	43 cms (17")
Maximum height	56 cms (22")
Diameter	32 cms (12")
Base length width	46 cms X 46 cms ( 18" x 18")

## OPTIONAL FEATURES AT EXTRA COST (ACCESSORIES)

• Five castors for stability.

#### STANDARD FEATURES

- Made of mild steel powder coated.
- Oven baked.
- Height adjustable with screw.
- Powder coated base with 4 nos. castors.
- Cushion top with rexine cover.

TECHNICAL SPECIFICATIONS	Approximate
Height maximum Height minimum	, ,
Overall width	4.4.5 (0.0.511)
Cushion sheet size: dia	36 cms (14")

## OPTIONAL FEATURES AT EXTRA COST (ACCESSORIES)

• Rubber shoes (4 nos.)

- Height adjustable with gas spring.
- Plastic non rusting base with five castors for stability.
- Revolving backrest revolves independently from seat.
- Cushioned seat and back.

TECHNICAL SPECIFICATIONS*	Approximate
Minimum height	51cms (20")
Maximum height	63 cms (25")
Length & width of heavy duty plastic base	
with non-rusting swivel castors	59cms (23")
Diameter of padded seat	36.20 cms (14.5")

### ATTENDANT CHAIR / BED (FOLDABLE)

P.C. No.: 1286

#### TRIPLE SEATING CHAIR

P. C. No.: 1290





#### STANDARD FEATURES

- Strong tubular construction.
- Padded seat back and sleep area.
- In built headrest.
- Easy folding arrangement.

TECHNICAL SPECIFICATION	ONS <sup>*</sup> Approximate
As a chair:	
Max height (back rest)	94 cms ( 37" )
Seat height	
Seat area	61 cms X 56 cms ( 24"X 22")
As a bed:	
Max height (including arm rest)	53.5 cms ( 21")
Sleep deck height	
Sleep deck area (L x B)	
Safe working load	9 .

- Strong rectangular base with black powder coated.
- $\bullet$  Tabular molded chrome plated leg & arm rest.
- Perforated powder coated seat & back rest.
- Every seat edge in chrome plated.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	180 cms (71")
Overall width	76 cms (30")
Overall height	78.5 cms (31")
Seat length	167.5 cms (66")
Seat width	45.5 cms (18")
Seat height	38 cms (15")



### **BLOOD DONOR / DIALYSIS COUCH (ELECTRIC)**

P.C. No.: 1277



#### STANDARD FEATURES

- Back rest, Knee rest and height adjustable on actuator.
- 3 actuators with wired hand set.
- 3 section of steel frame with padded mattress of rexine with one pillow.
- Telescopic I.V. pole, can be located at head ends.
- Urine bag holder provided on both sides.
- Height adjustable arm rest, with swing facility.
- Four castors 12.5 cms (5"), two castors with brakes and two without brakes.
- Epoxy, powder coated finished.
- Stainless steel 3 tired side try.
- Height adjustable fully S.S. foot rest.

#### TECHNICAL SPECIFICATIONS\*

Approximate

Overall length	186 cms (73")
Overall width with S.S. trays	
Overall width without S.S. trays	96.5 cms (38")
Minimum height	63 cms (25")
Maximum height	89 cms (35")
Seating area	51cms x 60 cms (20" x 23.5")
Lying area	186 cms x 60 cms (73" x 23.5")
Features from chair position to r	eclining & lying position
with wired hand remote	

### BLOOD DONOR / DIALYSIS COUCH (MANUAL)

P.C. No.: 1277M



- 3 manual crank on bearing for back knee and hi-lo.
- 3 section of steel frame with padded.
- 2 sectional mattress of rexine with pillow.
- Telescopic I.V. pole, can be located at head end.
- Urine bag holder provided on both side.
- Armrest with height adjustable, sliding type and with swing facility.
- 12.5 cms (5") dia plastic body castor.
- Epoxy powder coated finished.
- Height adjustable fully S.S. foot rest.

TECHNICAL SPECIFICATIONS	Approximate
Overall length Seating area Overall width	175 x 63 cms (69" x 25")
Minimum Height	59 cms (23")
Chair Position	
Back support area	63 x 76 cms (25" x 30")
Seat area	63 x 48 cms (25" x 19")