GYNAEC & PAEDIATRIC

1.



EXAMINATION COUCH GYNAEC EXAMINATION COUCH GYNAEC P.C. No.: 1242 Back raising without Gas Spring Back raising with Gas Spring

STANDARD FEATURES

- Made of M. S. with powder coating.
- Two lockable sliding doors on S. S. channel.
- Two pull drawers, one lockable.
- Heavy padded top with water resistant vinyl fabric with "U" cut.
- One stainless steel I.V. pole.
- One BP instrument tray.
- Adjustable and sliding lithomy calf rest.
- Stainless steel bowl.
- Built in step stool.

TECHNICAL SPECIFICATIONS^{*}


 Overall length
 198 cms (78")

 Overall width
 63.5 cms (25")

 Overall height
 84 cms (33")

Sleep deck area _____195.5 cms x 63.5 cms (77" x 25")

GYNAEC



STANDARD FEATURES

- 2 sectional top made from heavy 16 gauge pipe.
- 2 section padded Top with water Resistant vinyl fabric with "U" cut.
- Back section operated with gas spring OR manual height adjustable pipes.
- Adjustable and sliding lithomy calf rest on both side.
- Removable and sliding type S.S. bowl provided in leg side below padded seat.
- Legs made from 16 gauge M.S. tubular pipes.

TECHNICAL SPECIFICATIONS*

Over all length	
Over all width	63 cms (25")
Over all height	81 cms (32")
Sleep deck area	198 x 63 cms (78" x 25")



BIRTHING BED / LDRP BED (ELECTRIC)

P. C. No.: 1281E

Gynaec



STANDARD FEATURES

- 3 electrical function for back, hi-lo and T/RT which operated by wired hand set.
- 3 section of steel frame with 8 mm thick backelite heavy sheet top with telescopically sliding leg section.
- 3 sectional removable foam mattress provided on backelite top with seperate leg section.
- S.S. basinet Located at below fixed section and "C" cut given in fixed section mattress for basinet.
- Urine bag holder provided in both side.
- Central locking system for castors .
- S.S. Lithotomy crutches provided with swivel calf support.
- Below the deck sliding P.P. railing on both sides.
- P.P. head & foot board on both end.
- Battery back up provided for all 3 functions.
- Sliding bearing down hand grip provided for extra strength for patient.
- Epoxy powder coated finish.

TECHNICAL SPECIFICATIONS*

Overall Length with buffer	
Overall width with buffer	
Sleeping area	188 x 76 cm (74" x 30")
Maximum height	91 cm (36")
Minimum height	61 cm (24")

GYNAEC



STANDARD FEATURES

- Backrest on ratchet, head-low/up by screw mechanism.
- Middle section having "U"notch, with bowl.
- Below the deck telescopically sliding foot end.
- Railing on both sides.
- Lithotomy crutches with calf support .
- Bearing down grip handle.

OPTIONAL FEATURES AT EXTRA COST (ACCESSORIES)

- Step stool.
- SS telescopic I.V. pole.
- 35mm thick mattress.

TECHNICAL SPECIFICATIONS^{*}

Overall length	
Overall width	
Overall height	77 cms (30¼")
SS rectangular tube top frame	5 x 2.5 cms (2" x 1")
SS tubular bottom frame	3.2 cms (1¼")



BABY CRIB

P. C. No.: 1222 With Rod for Maternity Bed



STANDARD FEATURES

- Made of mild steel powder coated.
- Oven baked.
- MS stand for attachment with bed (powder coated).

TECHNICAL SPECIFICATIONS	Approximate
Overall length of crib	30" (76 cms)
Overall width of crib	
Over all height of crib	
Sleeping area of crib	14" to 28" (35.5 cms to 71
cms)	
Over all length with stand	
Overall width With stand	16" (41 cms)
Overall height of crib with stand	
from mattress frame	27" (68.5 cms)
Distance between crib and	
mattress frame of bed	12" (30.5cms)

PAEDIATRIC BED P. C. No.: 1224 For Children / Infants



STANDARD FEATURES

- Made of mild steel powder coated.
- Oven baked.

- With MS drop railings (powder coated).
- Four non rusting castors 12.5 cms (5"), two castors with brakes and two without brakes.

TECHNICAL SPECIFICATIONS^{*}

Overall length with railings	161 cms (63.5")
Overall width with railings	
Over all height with railings from	floor 106 cms (41.75")
Sleep deck area	
length & width	152 cms to 76 cms (60" x 30")
Height from floor to	
sleep deck with castors	61 cms (24")

PAEDIATRIC

BABY CRIB

P.C. No.: 1225 On S. S. Trolley



STANDARD FEATURES

- Made of stainless steel 304.
- Detachable moulded crib on trolley.
- Tilting reverse trendelenburg head side up.
- S. S. basket at bottom for storage.
- Four non rusting castors 12.5 cms (5"), two castors with brakes and two without brakes.
- 3 cms thick mattress with water resistant vinyl fabric.

TECHNICAL SPECIFICATIONS^{*}

Overall length	81 cms (32")
Overall width	51 cms (20")
Overall height with 12.5 cms (5") case	tors 108 cms (42.5")
Crib length	
Crib width	
Tilting reverse trendelenburg head	d side up 10°
Mattress length	68.5 cms (27")
Mattress width	37 cms (14 ½")
S. S. basket	45 cms x 48 cms (17.5" x19")
Height from floor level	25 cms (10")