



HOSPITAL BEDS (Electric)

PLATINUM CARE 4000

P.C. No.: 1200

Fully Electric ICU Bed with night lamp facility,
 Bed extender facility, PP Panels, Five electrically operated
 movements: Back, Hi-Low Knee, T/RT, Battery backup system,
 Angle indicator and base cover.



STANDARD FEATURES

- 4 actuators and ACP (Attendant Control Panel) at foot end.
- Non powered instant CPR with cable release on both sides.
- 4 sections of steel frame with mattress area of CRC strips of 18 gauge to facilitate cleaning and better ventilation as per FDA norms world wide.
- Anti squeeze back rest and knee movements.
- Easily removable moulded head and foot boards .
- Polymer moulded below the sleep deck retractable full length split railings.
- Four position I.V. pole mounts.
- One stainless steel telescopic I.V. pole with storage under sleep deck.
- Urine bag holder provided on both sides.
- Facility for lifting pole.
- Four non rusting castor with central locking.
- Oxygen cylinder cage.
- Corner roller buffers.
- Epoxy powder coated finished.

OPTIONAL FEATURES AT EXTRA COST (ACCESSORIES)

- Lifting pole.
- 4 sectional mattress.
- Top section in ABS.
- Base cover.
- X-ray permeable backrest with cassette holder.
- Perforated sheet top.

TECHNICAL SPECIFICATIONS*

Approximate

Overall bed length including buffers	228 cms (89.75")
Width with side rail up	104 cms (41")
Width with side rail down	99 cms (39")
Height of sleep deck (low position)	46 cms (18")
Height of sleep deck (high position)	71 cms (28")
Head section articulation	0 through 80°
Knee section articulation	0 through 35°
Trendelenburg and reverse trendelenburg	13°
Sleep deck area	204 x 90 cms (80" x 35.5")
4 section mattress size	201 x 89 cms (79" x 35")
Castor size approx	12.5 cms (5")
Under clearance when	
Sleep deck at height of 56 cms (22")	13 cms (5")
Safe working load	220 kgs (485 lbs.)

ULTIMA CARE 4000

P.C. No.: 1201

Fully Electric ICU Bed with CPR and Battery Backup System



STANDARD FEATURES

- 4 actuators with wired hand set and ACP (Attendant Control Panel) at foot end.
- Non powered instant CPR with cable release on both sides.
- 4 sections of steel frame with mattress area of CRC strips of 18 gauge to facilitate cleaning and better ventilation as per FDA norms world wide.
- Anti squeeze back rest and knee movements.
- Easily removable moulded head and foot boards .
- Polymer moulded below the sleep deck retractable full length split railings.
- Four position I.V. pole mounts.
- One stainless steel telescopic I.V. pole with storage under sleep deck.
- Urine bag holder provided on both sides.
- Facility for lifting pole.
- Four non rusting castor with central locking.
- Oxygen cylinder cage.
- Corner roller buffer.
- Epoxy powder coated finished.

TECHNICAL SPECIFICATIONS*

Approximate

Overall bed length including buffers	226 cms (89")
Width with side rail up	104 cms (41")
Width with side rail down	99cms (39")
Height of sleep deck (low position)	46 cms (18")
Height of sleep deck (high position)	71 cms (28")
Head section articulation	0 through 80°
Knee section articulation	0 through 35°
Trendelenburg and reverse trendelenburg	13°
Sleep deck area	204 x 90 cms (80" x 35.5")
4 section mattress size	201 x 89 cms (79" x 35")
Castor size approx	12.5 cms (5")
Under clearance when sleep deck at height of 56 cms (22")	13 cms (5")
Safe working load	220 kgs (485 lbs.)

OPTIONAL FEATURES AT EXTRA COST (ACCESSORIES)

- Lifting pole.
- 4 section mattress.
- X-Ray permeable backrest with cassette holder
- Base cover.
- Angle Indicator.
- Perforated sheet top.



UNIQUE CARE 2000

P.C. No.: 1205

Semi Electric ICU Bed with backrest and height adjustment with electric actuators with P. P. Panels



STANDARD FEATURES

- 2 actuators with wired hand set for back and hi-low which can be cranked for emergency with special manually operated crank handle.
- Non powered instant CPR with crank mechanism.
- 4 sections of steel frame with mattress area of CRC strips of 18 gauge to facilitate cleaning and better ventilation as per FDA norms world wide.
- Anti-squeeze backrest and knee movements.
- Moulded head and foot boards easily removable.
- Below the sleep deck collapsible railings.
- Four position I.V. pole mounts.
- One stainless steel telescopic I.V. pole with storage under sleep deck.
- Urine bag holder provided on both sides.
- Facility for lifting pole.
- Four non rusting castors 12.5 cms (5"), two castors with brakes.
- Corner roller buffers.
- Epoxy powder coated finished.

OPTIONAL FEATURES AT EXTRA COST (ACCESSORIES)

- Lifting pole.
- 4 section mattress.
- Oxygen cylinder cage.
- Moulded split railings, retractable below the sleep deck.
- Aluminium collapsible railing. (P.C. No. 1205)
- Split railing. (P.C. No.1205A)
- Base cover.
- Perforated sheet top.

UNIQUE CARE 2000

P.C. No.: 1205A

Semi Electric ICU Bed with ABS top, P.P. Panels and Height adjustment backrest with electric actuators



TECHNICAL SPECIFICATIONS*

Approximate

Overall bed length including buffers	226 cms (89")
Width with side rail up	104 cms (41")
Width with side rail down	104 cms (41")
Height of sleep deck (low position)	51 cms (20")
Height of sleep deck (high position)	71 cms (28")
Head section articulation	0 through 80°
Knee section articulation	0 through 35°
Trendelenburg and reverse trendelenburg	10°
Sleep deck area	204 x 90 cms (80" x 35.5")
4 section mattress size	201 x 89 cms (79" x 35")
Castor size	12.5 cms (5")
Under clearance when sleep deck at height of 56 cms (22")	13 cms (5")
Safe working load	220 kgs (485 lbs.)