



MULTIPURPOSE UTILITY PRODUCTS

INSTRUMENT TROLLEY

P.C. No.: 1254

Fully Stainless Steel



STANDARD FEATURES

- Made of stainless steel 304.
- Three shelves, top shelf with open front and railings on three sides.
- Middle and lower shelf, railings on all four sides.

TECHNICAL SPECIFICATIONS*

Approximate

Length	69 cms	(27")
Width	51 cms	(20")
Height with 12.5 cms (5") castors	91.5 cms	(36")
Length with 12.5 cms (5") castors push handle	77.5 cms	(30.5")
Distance between upper and middle shelf	25.5 cms	(10")
Distance between middle and lower shelf	25.5 cms	(10")

MEDICINE AND FOOD TROLLEY

P.C. No.: 1255

Fully Stainless Steel



STANDARD FEATURES

- Made of stainless steel 304.
- Two shelves, railings on all four sides.

TECHNICAL SPECIFICATIONS*

Approximate

Length	68 cms	(27")
Width	51 cms	(20")
Height with 12.5 cms (5") castors	91.5 cms	(36")
Length with 12.5 cms (5") castors and push handle	79 cms	(31")
Distance between upper and lower shelf	53 cms	(21")

DRESSING TROLLEY WITH BOWL AND BUCKET

P.C. No.: 1256

Fully Stainless Steel



STANDARD FEATURES

- Made of stainless steel 304.
- Two shelves, railings on all four sides.

TECHNICAL SPECIFICATIONS*

Approximate

Length	75 cms	(29")
Width	58 cms	(20")
Height with 12.5 cms (5") castors	91.5 cms	(36")
Length with 12.5 cms (5") castors, push handle and bowl & bucket	118 cms	(46.5")
Bowl diameter	30.5 cms	(12")
Bucket diameter	25.5 cms	(10")
Bucket height	20.5 cms	(9")
Distance between upper and lower shelf	53 cms	(21")

SOILED LINEN TROLLEY

P.C. No.: 1257S

Square



STANDARD FEATURES

- Made of stainless steel 304 quality.
- One heavy duty canvas bag.
- Four non-rusting, swivel castors 12.5 cms (5"), two castors with brakes and two without brakes.

TECHNICAL SPECIFICATIONS*

	Approximate
Length with handle	58.50 cms (23")
Width	47 cms (18.5")
Height with 12.5 cms (5") castors	103 cms (40.5")

SOILED LINEN TROLLEY

P.C. No.: 1257R

Round



STANDARD FEATURES

- Made of stainless steel 304 quality.
- One heavy canvas bag.
- Three non-rusting, swivel castors 12.5 cms (5"), one castor with brakes and two without brakes.

TECHNICAL SPECIFICATIONS*

	Approximate
Overall diameter of base	51cms (20")
Overall height	94cms (37")

MONITOR TROLLEY

P.C. No.: 1263

Deluxe



STANDARD FEATURES

- Made of mild steel powder coated, oven baked.
- 4 height adjustable shelves.
- One sliding pull drawer at bottom.
- One oxygen cylinder cage.
- One stainless steel I.V. pole, adjustable in height.
- Four castors 12.5 cms (5"), two castors with brake and two without brakes.
- 4 nos 5/15 amps electric points with one MCB.

TECHNICAL SPECIFICATIONS*

Approximate

Overall length	83 cms (33")
Length without cylinder cage	66cms (26")
Overall width	56 cms (22")
Overall height	155 cms (61")
Shelf size	51 cms x 45 cms (20"x18")
Usable drawer size	
Width	46 cms (18")
Depth	13 cms (5")
Pull	44 cms (17")

ECG TROLLEY

P.C. No.: 1287 With EPC frame work

P.C. No.: 1287S With S.S. frame work



With EPC frame work



With S.S. frame work

STANDARD FEATURES

- Made of mild steel powder coated.
- Oven backed.
- 5.0 cms heavy duty non-rusting swelving castors.
- I. V. pole with 2 hooks. (stainless steel)
- S.S. pull push handle or hand grip.

TECHNICAL SPECIFICATIONS*

Approximate

Overall length with S.S. handle	51 cms	(20")
Overall width with I.V. knob	44.5 cms	(17.5")
Overall height with S.S. handle	82.5 cms	(32.5")
Tray (shelf) size: length (upper shelf)	44.5 cms	(17.5")
Tray (shelf) size: width (upper shelf)	37 cms	(14")
Tray (shelf) size: middle & lower shelf (length)	44.5 cms	(17.5")
Tray (shelf) size: middle & lower shelf (width)	30.5 cms	(12")
Distance between lower & middle shelf	19 cms	(7.5")
Distance between upper & middle shelf	21.5 cms	(8.5")
Handle height from upper shelf	22.5 cms	(9")

TRAUMA CARE CRASH CART

P.C. No.: 1261

Old Design



STANDARD FEATURES

- Fully S.S. 304 grade.
- Two compartment with three drawers.
- Oxygen cylinder cage on one side.
- Wooden tray.
- Central lock for drawer.
- 6 plastic bins.
- Four castors 12.5 cms (5"), two castors with brake and two without brakes.

TECHNICAL SPECIFICATIONS*

	Approximate
Overall length	93 cms (36.5")
Overall width	56 cms (22")
Overall height	159 cms (62")

TRAUMA CARE CRASH CART

P.C. No.: 1261A

Deluxe with 4 Drawers



STANDARD FEATURES

- Fully S.S. 304 grade.
- Four compartment drawers.
- Oxygen cylinder cage on one side.
- Wooden tray.
- Central lock for drawer.
- 6 plastic bins.
- Four castors 12.5 cms (5"), two castors with brake and two without brakes.

TECHNICAL SPECIFICATIONS*

	Approximate
Overall length	60 cms (27.5")
Overall width with cylinder cage & handle	90 cms (35.5")
Overall height with top	162.5 cms (63.75")
Overall height without top	106.5 cms (42")
Overall box length	93 cms (36.5")
Overall box width	56 cms (22")
Overall box height	159 cms (62")

STAINLESS STEEL
OPERATION THEATRE TROLLEY

P.C. No.: 1262



STANDARD FEATURES

- Made from 304 stainless steel.
- Four side tray construction on top with 3 side railings and front open.
- One lower shelf with three side railing and front open.
- Four mounted on non-rusting castors 12.5 cms (5"), two castors with brake and two without brakes.

TECHNICAL SPECIFICATIONS*

Approximate

Overall length	91.5 cms (36")
Overall width	61 cms (24")
Overall height	96.5 cms (38")
Tray height with castors	86.5 cms (34")
Distance between upper tray and lower shelf	66 cms (26")