

A grayscale photograph of a hospital gurney, showing its metal frame, wheels, and adjustable bed surface. A red rectangular box is overlaid on the right side of the image, containing white text. The gurney is positioned diagonally, with its head towards the top left and its foot towards the bottom right. The background is a plain, light-colored wall.

PATIENT TRANSFER EQUIPMENT

HEIGHT ADJUSTABLE PATIENT STRETCHER

P.C. No.: 1258

On Trolley



STANDARD FEATURES

- Detachable two section stretcher top.
- Head rest can be raised manually on ratchet.
- Height adjustable on crank mechanism.
- Heavy padded ply top with water resistant vinyl fabric.
- Three built-in patient safety straps fixed on stretcher.
- One stainless steel I.V. pole which can be fitted on either side at head ends.
- Urine bag holder provided on both sides.
- Trolley provided with oxygen cylinder cage.
- 4 rotating round buffers to prevent wall or door damage.
- Four heavy duty castors 15.0 cms, two castors with brakes and two without brakes.
- Epoxy powder coated finished.
- Collapsible side railings.

TECHNICAL SPECIFICATIONS*

Approximate

Stretcher length with buffers	210cms (83")
Stretcher length	198cms (78")
Stretcher width	56cms (22")
Stretcher width with buffers	61 cms (24")
Trolley with stretcher minimum height	56 cms (22")
Trolley with stretcher maximum height	81 cms (32")
Trolley length	150 cms (59")
Trolley width with safety railings	68.5 cms (27")

HEIGHT ADJUSTABLE
PATIENT STRETCHER

P.C. No.: 1259

On Trolley with Longitudinal Tilting (Trendelenberg and Reverse Trendelenberg)



STANDARD FEATURES

- Detachable stretcher top.
- Head rest raised manually on ratchet.
- Height adjustable and longitudinal tilting (Trendelenburg and Reverse Trendelenburg) with 2 separate cranks on bearings.
- Trendelenburg and Reverse Trendelenburg $\pm 6^\circ$ approx.
- Heavy padded ply top with water resistant vinyl fabric.
- Three built in patient safety straps fixed on stretcher.
- One stainless steel I.V. pole which can be fitted on either side at head ends.
- Urine bag holder provided on both sides.
- Trolley provided with oxygen cylinder cage.
- 4 rotating round buffers to prevent wall or door damage.
- Four heavy duty castors 15.0 cms, two castors with brakes and two without brakes.
- Epoxy, powder coated finished.

TECHNICAL SPECIFICATIONS*


Approximate

Stretcher length with buffers	210 cms	(82.5")
Stretcher length	198 cms	(78")
Stretcher width with buffers	61 cms	(24")
Trolley with stretcher minimum height	56 cm	(22")
Trolley with stretcher maximum height	81 cms	(32")
Trolley length	150 cms	(59")
Trolley width with safety railings	68.5 cms	(27")
Stretcher width	56 cms	(22")

FOLDING STRETCHER

P.C. No.: 1260

Stainless Steel with Aluminium Legs



STANDARD FEATURES

- Made from stainless steel tube and aluminium legs.
- Legs to keep patient elevated from floor.
- Foldable along the width.
- Lockable brackets at the bottom during use. unlock for folding.
- Strong heavy duty vinyl fabric fitted with restraint straps for patient safety.
- Washable.

TECHNICAL SPECIFICATIONS*

Approximate

Stretcher length	223 cms	(88")
Stretcher width	53.5 cms	(21")

Foldable dimensions

Width	10 cms	(4")
Height	14 cms	(6")

PATIENT STRETCHER TROLLEY

P. C. No.: 1276

Fixed height, Aluminium top & M.S. frame
(Supplied with 25 mm mattress, Saline rod,
Oxygen cylinder cage)



STANDARD FEATURES

- Made of mild steel powder coated.
- Oven baked.
- With Aluminium top & MS frame with 25 mm thick mattress.
- With Oxygen cylinder cage.
- Four heavy duty castors (15.0 cms), two castors with brakes and two without brakes.
- Detachable stretcher on 6 legs on trolley.
- With fully SS saline rod.

TECHNICAL SPECIFICATIONS*

	Approximate
Overall length	203 cms (80")
Overall width	59.5 cms (23")
Overall height	82 cms (32.5")
Stretcher length & width	203 cms x 56 cms (80" x 22")
Sleep deck area	183cms x 56 cms (72" x 22")

OPTIONAL FEATURES AT EXTRA COST (ACCESSORIES)

- Swing away railings.

TRAUMA CARE EMERGENCY AND RECOVERY TROLLEY

P.C. No.: 1279 (Without Pump)

P.C. No.: 1279A (With Pump)



STANDARD FEATURES

- Detachable stretcher on 4 legs on trolley.
- Head rest adjustable on ratchet.
- Height adjustable and longitudinal tilt (Trendelenberg and reverse trendelenberg) on two cranks mechanism.
- X-Ray permeable stretcher top with removable padded vinyl cushion.
- X-Ray sliding tray along the length of the trolley.
- Stainless Steel (304 quality.) Safety railings on both sides, Swingable drop type.
- One stainless steel I.V. Pole can be located at two positions on head ends.
- Two protective stainless steel handles used also for lifting stretcher.
- Urine bag holder on both sides.
- Oxygen cylinder cage.
- Corner roller bumpers to prevent wall or door damage.
- Epoxy oven baked.

TECHNICAL SPECIFICATIONS*

	Approximate
Stretcher length with bumpers	210 cms (82.5")
Stretcher length without bumpers	198 cms (78")
Stretcher width with bumpers	61 cms (24")
Stretcher width without bumpers	61 cms (24")
Trolley length without railings	198 cms (78")
Trolley width with railings	66 cms (26")
Stretcher height minimum	73.5 cms (29")
Stretcher height maximum	99 cms (39")
Trendelenburg reverse trendelenburg	10° ± 2
Castors 4 Nos. 2 Nos. with brakes	15 cms(6")

O.T. TRANSFER TROLLEY

P.C. No.: 1280



STANDARD FEATURES

- Detachable stretcher on trolley with legs.
- Stretcher transfer with automatic safety locks.
- Removable foam padded vinyl cushion on stretcher.
- One stainless steel telescopic I.V. pole can be located at 4 positions on stretcher.
- Stretcher transfer wheels on bearings for ease in transfer.
- P.P. safety railings on stretcher.
- Protective and lifting S.S. handles on stretcher.
- Castors with central locking.
- Height adjustable trolleys.

TECHNICAL SPECIFICATIONS*

Approximate

Stretcher length with handle	206 cms	(81")
Stretcher width without railings	61 cms	(24")
Stretcher width with railings	85 cms	(33.5")
Travel trolley length	152 cms	(60")
Travel trolley width	61 cms	(24")
Trolley main frame width	48 cms	(19")
Trolley main frame length	168 cms	(66")
Total length of two trolley	355 cms	(140")
Stretcher height minimum	79 cms	(31")
Stretcher height maximum	107 cms	(42")