



ADJUSTABLE BEDSIDE TABLE

P.C. No.: 1244 Membrane top with crank mechanism
P.C. No.: 1246 Membrane top with gas spring mechanism



ADJUSTABLE BEDSIDE TABLE

P.C. No.: 1245 S.S. top with crank mechanism
P.C. No.: 1247 S.S. top with gas spring mechanism



STANDARD FEATURES

- Approximate size : L 75 cms (29.5") x W 40 cms (16")
- Height adjustable.
- Table top with adjustable angle for reading.
- M.S. tubular base frame mounted on non-rusting 5.0 cms swiveling twin castors.
- Pre-treated and powder coated.
- Supplied in KDC.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	75 cms (29.5")
Overall width	40 cms (16")
Minimum height	76 cms (30")
Maximum height	115 cms (45")

STANDARD FEATURES

- Approximate size : L 75 cms (29.5") x W 40 cms (16")
- Height adjustable.
- Table top with adjustable angle for reading.
- M.S. tubular base frame mounted on non-rusting 5.0 cms swiveling twin castors.
- Pre-treated and powder coated.
- Supplied in KDC.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	75 cms (29.5")
Overall width	40 cms (16")
Minimum height	76 cms (30")
Maximum height	115 cms (45")

ADJUSTABLE BEDSIDE TABLE NEW TYPE

P.C. No.: 1244M Membrane top with crank mechanism P.C. No.: 1244MG S.S. top with gas spring mechanism



- Approximate size : L 75 cms (29.5") x W 40 cms (16")
- Height adjustable.
- Table top with adjustable angle for reading.
- M.S. tubular base frame mounted on non-rusting 5.0 cms swiveling twin castors.
- Pre-treated and powder coated.
- Supplied in KDC.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	75 cms (29.5")
Overall width	40 cms (16")
Minimum height	76 cms (30")
Maximum height	115 cms (45")

MAYO TABLE

P.C. No.: 1248G

Height adjustable with gas spring



MAYO TABLE

P.C. No.: 1248C

Height adjustable with geared crank on bearing



FOLDING FOOD TABLE

P.C. No.: 1288

Fully stainless steel



STANDARD FEATURES

- Made of 304 stainless steel.
- \bullet 5.0 cms heavy duty non-rusting swelving castors.

TECHNICAL SPECIFICATIONS*	Approximate
for 1248 G	
Overall length	61 cms (24")
Overall width	40.5 cms (16")
Minimum height	73 cms (28.75")
Maximum height	103 cms (40.5")

STANDARD FEATURES

- Made of 304 stainless steel.
- 5.0 cms heavy duty non-rusting swelving castors.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	61 cms (24")
Overall width	40.5 cms (16")
Minimum height	76 cms (30")
Maximum height	114 cms (45")

- Made of stainless steel (304).
- 4 nos. rubber shoes.

TECHNICAL SPECIFICATIONS*	Approx	kimate
Overall length	56 cms	(22")
Overall width	42 cms	(16.5")
Overall height with rubber shoes	28 cms	(11")
Foldable height	11.5 cms	(4.5")



ECONOMY BEDSIDE LOCKER

P.C. No.: 1250

With Membrane Top



ECONOMY BEDSIDE LOCKER

P.C. No.: 1251

With Stainless Steel Top



STANDARD FEATURES

- Made of marine ply laminated top.
- One magnetic catch door.
- One open shelf.
- Optional with or without castors.
- Epoxy powder coated finished.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	40.5 cms (16")
Overall width	37 cms (14.5")
Overall height	81 cms (32")

- Made of stainless steel top.
- One magnetic catch door.
- One open shelf.
- Optional with or without castors.
- Epoxy powder coated finished.

TECHNICAL SPECIFICATIONS	Approximate
Overall length	37.5 cms (14.75")
Overall width	33 cms (13")
Overall height	81 cms (32")
Laminated top net. weight	12.000 kgs
Stainless steel top net. weight	11.000 kgs

DELUXE BEDSIDE LOCKER

P.C. No.: 1252 With Membrane Top

P.C. No.: 1252M With Modified Membrane Top

DELUXE BEDSIDE LOCKER

P.C. No.: 1253

With Stainless Steel Top







With Modified Membrane Top



STANDARD FEATURES

- Locker made of machine pressed CRC steel sheets enclosed on three sides.
- One pull drawer on top. One open shelf.
- One magnetic catch door with inbuilt storage inside at bottom.
- Mounted on non-rusting 5.0 cms swiveling twin castors, two castors with brakes and two without brakes.
- Pretreated and powder coated.
- Buffer at front and back.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	41 cms (16")
Overall width	44 cms (17.5")
Overall height	86 cms (34")

- Locker made of machine pressed CRC and steel sheets enclosed on three sides.
- One pull drawer on top. One open shelf.
- One magnetic catch door with inbuilt storage inside at bottom.
- Mounted on non-rusting 5.0 cms swiveling twin castors, two with brakes.
- Pretreated and powder coated.
- Buffer at front and back.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	46 cms (18")
Overall width	43 cms (17")
Overall height	86 cms (34")



SINGLE STEP STOOL

P.C. No.: 1267

With Aluminium Sheet Top



STANDARD FEATURES

- Made of mild steel tube and sheets and covered with anti skid Aluminium sheet.
- Epoxy powder coated finished.
- 4 rubber shoes.

TECHNICAL SPECIFICATIONS*	Approximate
Step length	51 cms (20")
Step width	31 cms (12")
Step height	25.5 cms (10")

DOUBLE STEP STOOL

P.C. No.: 1268

With Aluminium Sheet Top



STANDARD FEATURES

- Made of mild steel tube and sheets and covered with anti skid Aluminium sheet.
- Epoxy powder coated finished.
- 4 rubber shoes.

TECHNICAL SPECIFICATIONS"	Approximate
Step length	51 cms (20")
Step width	23 cms (9")
First step height from floor	23 cms (9")
Second step height from floor	45 cms (18")

PATIENT STOOL

P.C. No.: 1269

With SS Sheet Top



- Made of mild steel tubes with SS Top.
- 4 Rubber shoes.
- Epoxy powder coated finished.

TECHNICAL SPECIFICATION	NS* Approximate
Square top (L x W):	30.5 cms X 30.5 cms (12" x 12"
) Square base (L x W): 38	3.5 cms X 38.5 cms (18.5" x 18.5")
Height of stool:	46 cms (18")

STAINLESS STEEL KICK BUCKET

P.C. No.: 1273



STAINLESS STEEL BOWL STAND

P.C. No.: 1274



STAINLESS STEEL BOWL STAND

P.C. No.: 1274A



STANDARD FEATURES

- Made of stainless steel.
- Three non rusting plastic body castors 12.5 cms (5"), two castors with brakes and one without brake.

TECHNICAL SPECIFICATIONS*	Approximate
Bowl diameter	38 cms (15")
Height	25 cms (10")

STANDARD FEATURES

- Made of stainless steel (304).
- Leg fitted with rubber shoes.
- 2 tier (bowls).

TECHNICAL SPECIFICATIONS*	Approximate
Overall height	91.5 cms (36")
Overall diameter	40 cms (16")
Bowl diameter	38 cms (15")
Distance between 2 bowls	36 cms (14")

OPTIONAL FEATURES AT EXTRA COST (ACCESSORIES)

• Three non rusting plastic body castors 12.5 cms (5"), two castors with brakes and one without brake.

STANDARD FEATURES

- Made of stainless steel (304).
- Leg fitted with rubber shoes.
- 1 tier (bowl).

TECHNICAL SPECIFICATIONS*	Approximate
Overall height	91.5 cms (36")
Overall diameter	40 cms (16")
Bowl diameter	38 cms (15")

OPTIONAL FEATURES AT EXTRA COST (ACCESSORIES)

• Three non rusting plastic body castors 12.5 cms (5"), two castors with brakes and one without brake.



I. V. STAND

P. C. No.: 1265

Regular (Epoxy, powder coated finished)



STANDARD FEATURES

- Five heavy duty plastic base castors for stability. (Non rusting)
- 4 hooks provided for I. V. fluids.
- Made of MS powder coating

TECHNICAL SPECIFICATIONS*	Approx	imate
Base length with castors	60 cms	(24")
Hook length	38 cms	(16")
Hook height minimum	136 cms ((53.5")
Hook height maximum	216 cms	(85")

I. V. STAND

P. C. No.: 1266

Fully Stainless Steel



- Five heavy duty S.S. base castors for stability. (Non rusting)
- 4 hooks provided for I. V. fluids.
- Made of stainless steel 304.

TECHNICAL SPECIFICATIONS	Approx	imate
Base length with castors	60 cms	(24")
Hook length	38 cms	(16")
Hook height minimum	136 cms (53.5")
Hook height maximum	216 cms	(85")

LEAD APRON STAND

P. C. No.: 1289



STANDARD FEATURES

- Made of mild steel powder coated.
- Oven baked.
- Five hooks with weight bearing capacity 6.5 kgs. each.
- Four heavy duty swelving castors 12.5 cms (5").

TECHNICAL SPECIFICATIONS*

Approximate

Overall length	109 cms (43")
Overall width	91 cms (36")
Overall height with castors	129.5 cms (51")
5 Hooks with length	57 cms (22.5")



WARDCARE OVER BED TABLE

P.C. No.: 1275

WARDCARE SINGLE OXYGEN CYLINDER TROLLEY P.C. No.: 1278





STANDARD FEATURES

- Made of mild steel.
- Laminated board on top.
- Epoxy powder coated finished.

TECHNICAL SPECIFICATIONS	Approximate
Overall length	113 cms (44.5")
Overall width	79 cms (31")
Laminated top	76 cms x 46 cms (30"x18")
Fixed height	97 cms (38")
Castors 4 nos	5 cms (2")

- Made of mild steel powder coated.
- Oven baked.

TECHNICAL SPECIFICATIONS*	Approximate
Overall length	35 cms (14")
Overall width	35 cms (14")
Overall height	128 cms (51")
Castors 2 nos	12.5 cms (5")

WARDCARE SCREEN

P.C. No.: 1264

With Fabric and Six Castors



STANDARD FEATURES

- Made of M. S. tubes with 3 folds.
- Folds on 6 castors.
- Fabric fixed on stainless steel rods.
- Easy to cleaning removable Fabric.
- Epoxy powder coated finished.

TECHNICAL SPECIFICATIONS* Approximate Screen height 180 cms (70")

243 cms (96")

121 cms (48")

Open screen length

Folded screen length